

**PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT  
FOR FUTURE USE IN CASE YOU HAVE A PROBLEM.**

## **BAUSCH & LOMB** Sensitive Eyes® Drops

**Bausch & Lomb® Sensitive Eyes Drops**, made only by Bausch & Lomb, is a gentle, effective thimerosal-free preserved solution that is recommended for rewetting your lenses during wear.

A virtually reaction-free formula that soothes eyes and moistens lenses to minimize dryness during wear. Ideal for rewetting lenses upon awakening for extended wear lens wearers.

**CONTENTS:** Bausch & Lomb Sensitive Eyes Drops is a sterile, buffered, aqueous solution that contains boric acid, sodium borate and sodium chloride; preserved with sorbic acid (0.1%) and edetate disodium (0.025%).

**ACTIONS:** Bausch & Lomb Sensitive Eyes Drops soothes eyes and moistens lenses to minimize dryness during wear.

**INDICATIONS (Uses):** Bausch & Lomb Sensitive Eyes Drops is indicated for use to rewet your lenses during wear. Also indicated for rewetting your lenses after awakening for extended wear lens wearers.

**CONTRAINDICATIONS (Reasons not to use):** If you are allergic to any ingredients in Bausch & Lomb Sensitive Eyes Drops, do not use this product.

**WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.** It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision.

Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight.

Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

Studies have also shown that smokers have a higher incidence of adverse reactions.

If you experience eye discomfort, excessive tearing, vision changes, or redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner.

It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

To avoid contamination, do not touch tip of container to any surface. Replace cap after using.

### **PRECAUTIONS:**

- Keep the bottle tightly closed when not in use.
- Store at room temperature.
- Use before expiration date marked on the bottle and carton.

**ADVERSE REACTIONS (Problems and what to do):** The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

- **IMMEDIATELY REMOVE YOUR LENSES.**
- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, **DO NOT** put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.
- If the problem continues, **IMMEDIATELY** remove the lens and consult your eye care practitioner.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

### **GOOD LENS CARE PRACTICES:**

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse lens case with fresh rinsing solution, and allow to air dry.

### **DIRECTIONS:**

- During wear, put two to three drops directly onto each lens whenever needed. Blink several times and wait a few moments.
- If the lens still does not feel comfortable, add another drop.
- If discomfort persists, **IMMEDIATELY** remove your lenses and contact your eye care practitioner.

**HOW SUPPLIED:** Bausch & Lomb Sensitive Eyes Drops is available for purchase in sterile 0.5 Fl. Oz. (15ml) and 1 Fl. Oz. (30ml) plastic bottles. Bottles and cartons are marked with a lot number and expiration date.

**PRODUCT SATISFACTION GUARANTEE**  
*If you are not completely satisfied with this product, call the toll-free number on this insert.*

**Toll Free Product Information**  
**Call: 1-800-553-5340**

**Bausch & Lomb** and **Sensitive Eyes** are registered trademarks of Bausch & Lomb Incorporated.

© Bausch & Lomb Incorporated

Manufactured by:  
Bausch & Lomb  
Rochester, NY 14692-0450  
Made in USA